

# Contact Lens MULTIPURPOSE SOLUTION

Specialized company in  
contact lens care solutions



Contact us!

Taewon Enterprise  
www.taewonenterprise.com

## INDEX

---

Eye+ Multi Solution  
Eye Premium Solution  
Pink Angel Multi Solution  
Saline Solution Eyeyou  
Clean Saline Solution Lloyds  
LK Cleaner Solution DL Plus  
Solution  
Rewetting Solution  
OFresh Resense Solution  
OFresh Solution  
Purecle Solution  
Global

Contact lens care all at once!

# EYE+ MULTI

## Multi-Purpose Solution

### HARD SOFT CONTACT LENS MULTIPURPOSE SOLUTION

- Excellent lubrication with hyaluronic acid
- Keeping Moisturized throughout the day from morning to afternoon
- Cleaning, disinfection, rinsing, preservation, protein removal, lens stabilization
- All-in-One lens solution
- Silicone hydrogel, continuous, daily, color, RGP lenses



EYE+ MULTI  
360ml / 150ml 

## Eye+ Multi Solution

All-in-One multi lens solution



### Efficacy

What this product does: It cleans and removes film, protein, and other additives or residues adhering to contact lenses.

With daily use of this product, protein can be effectively removed. It also kills harmful microorganisms on the surface of the lens. Also, it can be used for rinsing, storing, or preventing the lenses from drying out before wearing them.

It can be used for daily care of contact lenses(hydrophilic). It is also used for cleaning, deproteinizing, rinsing, chemical sterilization(not heat sterilization) and preservation.

### Usage

- Step 1: Drop 3 drops of this solution on both surfaces of the lens and rub for 20 seconds.
- Step 2: Rinse the surface of the lens thoroughly with this solution to remove any residue remaining on the surface of the lens.
- Step 3: Put the cleaned lenses in the lens case and soak them in this clean solution for at least 4 hours. These lenses can be worn at any time. There is no need to separately rinse again with saline. If any residue is left on the lens, they can be used after rinsing again with this solution.

The lens case cap should be closed if the lenses are not to be put on immediately.

Since saline does not have a sterilizing effect, do not mix this solution with saline.

### Protein removal from contact lenses

When cleaning and disinfecting lenses daily with this solution, proteins are also removed, so there is no need to use a separate protein remover or daily cleaning agent.

Ending lens care for only  
**Multi Solution**



# All-in-One

Both for Hard and Soft Lens

# EYE PREMIUM SOLUTION

Moist lenses, comfortable to wear  
No rubbing comfort



### Total Care System


- Cleaning
  - Rinsing
  - Disinfecting
  - Storing
  - Protein Removing
  - Moisturizing
- **Efficacy**  
Cleaning, disinfection, rinsing, Preservation, wetting, protein removal (soft, hard lens)

Eye Premium Solution 360ml / 150ml

PINKANGEL Lens Care Multi Solution

# Pink Angel Multi Solution

4 functions of lens care solution at once!



- **Product name**  
Pink Angel Multi Solution  
(Polyhexamethylene biguani hydrochloride solution 20%)
- **Efficacy**  
Cleaning, rinsing, disinfecting, protein removal and preservation, and wetting of soft contact lenses
- **Usage**  
4-hour usage (no rubbing and no rinsing)  
Put the lens in the lens case, fill it with this solution, close the lid, and leave it there for at least 4 hours.  
**5 minute usage**  
1 Drop 3 or more drops of this solution on each surface of the lens and rub gently for 10 seconds.  
2 Rinse the surface of the lens with this solution.  
3 Put the lens in the lens case, fill it with this solution, close the lid, and leave it there for more than 5 minutes. Lenses treated in this way can be worn directly on the eyes without rinsing separately. When using this solution for lens preservation, it can be stored for up to 30 days after putting the lens in the lens case, filling it with this solution, and closing the lid.

Pink Angel Multi Solution 360ml / 150ml

Saline Solution

# Saline Solution

**Always clean and comfortable Saline solution**  
Hard Lens • Dream Lens • Soft Lens  
Available for all lenses



Saline Solution 500ml

- Efficacy**  
Cleaning, preservation and chemical disinfection
- Usage**  
Use an appropriate amount as needed.

Saline solution for Contact lenses



# eyeyou

Contactlens Saline Solution  
**Eye you solution (Sodium chloride)**  
Celan and comfortable with biotechnology for eye  
EYEU CONTACTLENS SALINE SOLUTION

## EYEU CLEAN SALINE SOLUTION

Always clean and comfortable saline solution

- Used for cleaning and rinsing contact lenses
- Cleaner with automation of the entire process



Eyeyou Saline Solution 500ml

- Efficacy**  
Rinsing after general rinsing, thermal and chemical disinfection of contact lenses
- Usage**
  - General Rinsing: Always clean the lenses with lens cleaner and rinse with this solution before wearing or disinfecting the lenses.
  - Rinsing after thermal and chemical disinfection: After thermal and chemical disinfection, rinse thoroughly with this solution.

## RGP Contact Lens Cleaner

## Lloyds LK Cleaner Solution

### Cleaning Solution for RGP (Oxygen Permeable Hard Lens)

RGP lens cleaning more thoroughly

**Lloyds LK Cleaner Solution, only for hard lenses**

- **DAILY Cleaner Solution**  
For Rigid Gas Permeable and Hard Contact Lenses
- **Product Name**  
Lloyds LK Cleaner Solution(Miranol 2MCA Modified)
- **Efficacy**  
Cleaning of hard and oxygen-permeable hard lenses
- **Volume**  
50ml , 30ml
- **Usage**  
Add 3-4 drops to both sides of the lens and rub for about 20 seconds, then rinse both sides thoroughly with the rinsing solution.



Lloyds LK Solution 50ml, 30ml

## RGP Contact Lens Multipurpose

## DL Plus Solution

### Disinfection, cleaning, rinsing, preservation, and wetting of RGP lenses

One-stop RGP contact lens care

**Moist and comfortable DL Plus Solution**

- **RGP Contact Lens Multi Solution**  
It can disinfect, clean, rinse, preserve and even wet the lenses in one step. It is a colorless and transparent solution of the main component of polyhexamethylenebiguanide hydrochloride solution. It can be used with oxygen permeable(RGP, hard) lenses as well as dream lenses. In addition, it contains HPMC(Hydroxypropyl Methylcellulose), which maintains moisture for a long time with hyaluronic acid and greatly strengthens the moisture, so that the eyes are more comfortable and moist when using this solution. Use the DL Plus Solution, which contains all the functions of RGP lens care, for easy and convenient in your busy life.



DL Plus solution 150ml

Rewetting Solution

# Rewetting Solution for contact lenses

It alleviates the feeling of foreign substances when wearing contact lenses with rewetting solution



15ml Rewetting solution

## Eye Premium Solution

Relieving eye strain and comfortable to wear Drop 2 or 3 drops on the lens

- It is a rewetting solution that moisturizes the lenses for a long time by strengthening the moisturizing power, increases the lubrication of the lenses, and reduces the feeling of foreignness in the eyes.



15ml Rewetting solution

## Pink Angel Multi Solution

Rewetting Solution for Contact Lenses  
Moisturizing and comfortable

- Moisturizing throughout the day
- Moisturizing dry eyes and improving discomfort
- Soft lenses only



15ml Rewetting solution

## EYE PLUS MULTI

Rewetting Solution for Contact Lenses  
Drop 2, 3 drops on the lens!!  
Eye Plus Multi Solution

- Excellent lubrication with hyaluronic acid
- Continuously maintain the wetting of the lens surface
- Both for contact hard lenses and soft lenses

Domestic lens shop franchise OEM

Domestic lens shop franchise OEM



# OFresh Resense Solution

# OFresh Solution

**A total multi-solution for daily moist lenses**  
Moisturizing, comfortable and clean

**Lens moisturizer for daily moist lenses**  
Coat the moisture with a natural moisturizing factor



**Efficacy**  
Cleaning, disinfecting, rinsing, preserving, wetting and protein removal of contact lenses



If you feel dryness even while wearing the lenses, remove the lenses and wash them, then use a lens moisturizer.

OFresh Resense Solution 250ml

OFresh Solution 30ml



Domestic lens shop franchise OEM

# Purecle Solution

**Lens moisturizer for daily moist lenses**  
**Moisturizing, comfortable and clean**



- Efficacy**  
 Cleaning, disinfection, rinse preservation, wetting, protein removal of soft contact lenses and silicone lenses

Purecle Solution 150ml



## OFresh Resense Solution

Contact lenses

**Washing • Disinfecting • Rinsing • Preserving • Wetting and Protein Removal**

- Moisturize the lens for a long time by strengthening the moisturizing power!
- Increases lens lubrication and reduces feeling of foreign substances
- Wearing lenses as comfortable as my own eyes



## OFresh Solution

Contains hyaluronic acid to keep moisture for a long time

**Comfortable feeling without feeling of foreign substances**

- Moisturize the lens for a long time by strengthening the moisturizing power!
- Increases lens lubrication and reduces feeling of foreign substances
- Wearing lenses as comfortable as my own eyes



## Purecle Solution

Soft and silicone lenses

**Washing • disinfecting • rinsing • preserving • wetting • protein removal**

- Long-lasting moisture
- Safe lens care for the cornea
- Effective sterilization function
- Easy protein removal



## A global company that exports goods to the world

Developed and produced various lens care solutions of excellent quality with excellent price competitiveness, and exports them to the global market as well as the domestic market.

**(USA, Mexico, Vietnam, Kuwait, Iran, Paraguay, United Arab Emirates, Indonesia, and Saudi Arabia, etc.)**

In addition, we are supplying products under contract with leading lens specialty stores in Korea, such as O-Lens, Lens Town, and Lens Diva.



Certificate of  
Pharmaceutical  
(Quasi drugs)  
Manufacturing  
License



CE Certificate  
(MTIC)



ISO13485(IGC)



Certificate of  
Excellent  
Technology  
Evaluation  
Company



Certificate of R&D  
Department



Medical Device  
Manufacturing  
License